

ORTHO
SYSTEMS

2000
2026

Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Stehbrett Standing Platform



Features

Sizes available for the **smallest to large** patients

Unrivalled **20-month guarantee**

Engraved **serial number** for traceability

Clearly illustrated instructions for use

+ Easily adapted to the **next size Stehbrett** at any time

+ **Silicone bumper** to avoid damages and protect against injury

+ **Wheel locks** keep podium stationary

Standard

Fully assembled for clear cost estimates and uncomplicated ordering

+ **Pyramid adapters** for precise adjustment, up to 10°


+ **Easy removal** of the body shell

+ **Height adjustable** to accommodate growth

+ **Secure 3-point attachment** of the body shell at the feet and pelvic area

+ **Universal connecting plate** is interchangeable with the Leicht-Walker and Stehbrett

Upon request an **extra-wide connecting plate** can be supplied for more abduction





Article No.	Height Classification	Base (l x w)	Connecting Plate	Heel to Heel (middle)
Stehbrett-klein	65 cm to 100 cm	40 cm x 40 cm	25 cm x 16 cm	23 cm to 8 cm
Stehbrett-mittel	80 cm to 130 cm	45 cm x 45 cm	32 cm x 19 cm	30 cm to 12 cm
Stehbrett-groß	120 cm to 160 cm	50 cm x 50 cm	39 cm x 22 cm	36 cm to 18 cm
Stehbrett-groß-breit	120 cm to 160 cm	50 cm x 70 cm	39 cm x 22 cm	36 cm to 18 cm

With anti-tip protection to prevent toppling backwards



Article No. Combi-product	Height Classification	Base (l x w)	Connecting Plate	Heel to Heel (middle)
Stehbrett-KS-klein	65 cm to 100 cm	40 cm x 40 cm	25 cm x 16 cm	23 cm to 8 cm
Stehbrett-KS-mittel	80 cm to 130 cm	45 cm x 45 cm	32 cm x 19 cm	30 cm to 12 cm
Stehbrett-KS-groß	120 cm to 160 cm	50 cm x 50 cm	39 cm x 22 cm	36 cm to 18 cm
Stehbrett-KS-groß-breit	120 cm to 160 cm	50 cm x 70 cm	39 cm x 22 cm	36 cm to 18 cm

With anti-tip protection to prevent toppling backwards



Article No. Combi-product	Height Classification	Base (l x w)	Connecting Plate	Heel to Heel (middle)
Stehbrett-HG-klein	65 cm to 100 cm	40 cm x 40 cm	25 cm x 16 cm	23 cm to 8 cm
Stehbrett-HG-mittel	80 cm to 130 cm	45 cm x 45 cm	32 cm x 19 cm	30 cm to 12 cm
Stehbrett-HG-groß	120 cm to 160 cm	50 cm x 50 cm	39 cm x 22 cm	36 cm to 18 cm
Stehbrett-HG-groß-breit	120 cm to 160 cm	50 cm x 70 cm	39 cm x 22 cm	36 cm to 18 cm

With hand grips and anti-tip protection




Article No. Combi-product	Height Classification	Base (l x w)	Connecting Plate	Heel to Heel (middle)
Stehbrett-HG-KS-klein	65 cm to 100 cm	40 cm x 40 cm	25 cm x 16 cm	23 cm to 8 cm
Stehbrett-HG-KS-mittel	80 cm to 130 cm	45 cm x 45 cm	32 cm x 19 cm	30 cm to 12 cm
Stehbrett-HG-KS-groß	120 cm to 160 cm	50 cm x 50 cm	39 cm x 22 cm	36 cm to 18 cm
Stehbrett-HG-KS-groß-breit	120 cm to 160 cm	50 cm x 70 cm	39 cm x 22 cm	36 cm to 18 cm


Body Shell with Connecting Plate

Custom body shell used in conjunction with the Stehbrett supports primary deformities and selectively counters any possible secondary deformities.


Only a precise footbed and proper support of the ankle joint enable correct positioning of the knee and hip joints. the knee and hip joints.




Minimum length from connecting plate to 4-point grip:
30 cm (small)
36 cm (medium and large)




Growable




Stability



Protective Bumper




3-Point Fixation




Removable



Postural Control Exercises



HG = Hand Grips



KS = Anti-tip Protection

www.ortho-systems-inmotion.com