

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Foot-Walker Walking Podium




Features

Sizes available for the **smallest to large** patients

Robust aluminum construction

Standard

Fully assembled for clear cost estimates and uncomplicated ordering




Upon request an **extra-wide connecting plate** can be supplied for more abduction

Unrivalled **20-month guarantee**


Engraved **serial number** for traceability

Clearly illustrated instructions for use




Article No.	Height Classification	Connecting Plate	Heel to Heel (middle)
Foot-Walker-klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Foot-Walker-mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Foot-Walker-groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

With anti-tip protection to prevent toppling backwards




Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Foot-Walker- KS -klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Foot-Walker- KS -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Foot-Walker- KS -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

With hand grips



Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Foot-Walker- HG -klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Foot-Walker- HG -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Foot-Walker- HG -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

With hand grips and anti-tip protection




Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Foot-Walker- HG-KS -klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Foot-Walker- HG-KS -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Foot-Walker- HG-KS -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

Body Shell with Connecting Plate

Custom body shell used in conjunction with the Stehbrett supports primary deformities and selectively counters any possible secondary deformities.


Only a precise footbed and proper support of the ankle joint enable correct positioning of the knee and hip joints.




Minimum length from connecting plate to 4-point grip:
30 cm (small)
36 cm (medium and large)




Noise Dampening




Growable



Vertical Axis




Step Width



Step Length



3-Point Fixation



Removable



Reduced Weight



HG = Hand Grips



KS = Anti-tip Protection



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com