

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Leicht-Walker Walking Podium




Features

Sizes available for the **smallest to large** patients

Lightweight **carbon construction**
30% weight reduction

- + **Vertical alignment adjustment in four directions** without disassembly
- + **Parallel adjustment of step width** without disassembly
- + **Individual adjustment of step length** without disassembly
- + **Noise dampening movement stops**



Standard

Fully assembled for clear cost estimates and uncomplicated ordering

- + **Pyramid adapters** for precise adjustment, up to 10°
- + **Easy removal** of the body shell
- + **Height adjustable** to accommodate growth
- + **Secure 3-point attachment** of the body shell at the feet and pelvic area

Upon request an **extra-wide connecting plate** can be supplied for more abduction (aluminium)

- + **Universal connecting plate** is interchangeable with the Foot-Walker and Stehbrett

Unrivalled **20-month guarantee**


Engraved **serial number** for traceability

Clearly illustrated instructions for use



Article No.	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker-klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker-mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker-groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

With anti-tip protection to prevent toppling backwards



Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker- KS -klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker- KS -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker- KS -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

With hand grips



Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker- HG -klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker- HG -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker- HG -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

With hand grips and anti-tip protection




Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker- HG-KS -klein	65 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker- HG-KS -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker- HG-KS -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm











Body Shell with Connecting Plate


Custom body shell used in conjunction with the Stehbrett supports primary deformities and selectively counters any possible secondary deformities.

Only a precise footbed and proper support of the ankle joint enable correct positioning of the knee and hip joints.



Minimum length from connecting plate to 4-point grip:
30 cm (small)
36 cm (medium and large)





Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com