

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026



















Functional Overview Locked Knee Joints

All knee joints are shipped fully assembled. The joints are free of play, self-lubricating, dirt-resistant and available in at least 6 sizes per model.


Modular extension assist options can be added or changed at any time. These accessories enhance the function of the knee joint.



Orthosis	Locked	Extension	Flexion	Posterior Offset	Modular	Force Options	No Grinding Necessary	Locking Mechanism	Water Resistant	
										Varo-KA-0 Pages 108 -109
										Varo-PS-0 Pages 114 -115
										Konso-5 Pages 120 -121
										Konso-0 Pages 126 -127
										Klassik-Konso Pages 132 -133
										Blocker Pages 134 -135
										Klassik-Konso with Slide Pages 136 -137


Orthosis	Locked	Extension	Flexion	Posterior Offset	Modular	Force Options	No Grinding Necessary	Locking Mechanism	Water Resistant	
										 Varo-KA-0 with Dynamic-0-S Pages 110 - 111
										 Varo-KA-0 with Dynamic-0-N Pages 112 - 113
										 Varo-PS-0 with Dynamic-0-S Pages 116 - 117
										 Varo-PS-0 with Dynamic-0-N Pages 118 - 119
										 Konso-5 with Move-5-S Pages 122 - 123
										 Konso-5 with Move-5-N Pages 124 - 125
										 Konso-0 with Move-0-S Pages 128 - 129
										 Konso-0 with Move-0-N Pages 130 - 131

Locked Knee Joints with Extension Assist



Knee Joint Type	Varo-KA-0	Varo-KA-0	Varo-KA-0	Varo-PS-0	Varo-PS-0	Varo-PS-0	
Application	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	
Modular Add-on		Dynamic-0-S	Dynamic-0-N		Dynamic-0-S	Dynamic-0-N	
Force Adjustment		3 Positions	7 Positions		3 Positions	7 Positions	
Max. Flexion Angle	130°	115°	95 to 125°	130°	115°	95 to 125°	
Assitive Force		up to max. force	normal force		up to max. force	normal force	
Cylinder Detachable		possible			possible		
Locking Mechanism	Cable Release	Cable Release	Cable Release	Pre-select	Pre-select	Pre-select	
Unlock	Satellite	Satellite	Satellite	Push button	Push button	Push button	
Lock	Satellite	Satellite	Satellite	Push button	Push button	Push button	
Unlocked (maintained)	Satellite (unlocked)	Satellite (unlocked)	Satellite (unlocked)	Push button (unlocked)	Push button (unlocked)	Push button (unlocked)	
Donning Aid		yes	yes		yes	yes	Donning Aid Rudi
	Pages 108-109	Pages 110-111	Pages 112-113	Pages 114-115	Pages 116-117	Pages 118-119	Pages 146-147

► See pages 82–83
for extension assist
applications



Knee Joint Type	Konso-5	Konso-5	Konso-5	Konso-0	Konso-0	Konso-0	Klassik-Konso
Application	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO	KAFO/HKAFO/RGO
Modular Add-on		Move-5-S	Move-5-N		Move-0-S	Move-0-N	
Force Adjustment		3 Positions	7 Positions		3 Positions	7 Positions	
Max. Flexion Angle	130°	115°	95 to 125°	130°	115°	95 to 125°	130°
Assitive Force		up to max. force	normal force		up to max. force	normal force	
Cylinder Detachable		possible			possible		
Locking Mechanism	Drop lock	Drop lock	Drop lock	Drop lock	Drop lock	Drop lock	Drop lock
Unlock	Lift drop lock	Lift drop lock	Lift drop lock	Lift drop lock	Lift drop lock	Lift drop lock	Lift drop lock
Lock	Auto-locking	Auto-locking	Auto-locking	Auto-locking	Auto-locking	Auto-locking	Auto-locking
Unlocked (maintained)	Spring plunger	Spring plunger	Spring plunger	Spring plunger	Spring plunger	Spring plunger	(Spring plunger)
Donning Aid Rudi		yes	yes	yes	yes	yes	
	Pages 120-121	Pages 122-123	Pages 124-125	Pages 126-127	Pages 128-129	Pages 130-131	Pages 132-133

www.ortho-systems-inmotion.com