

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics


Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Addition Separator-S

Simple Operation in 3 Steps






Partially integrated into the upright bar, this quick-release system separates a **KAFO** into **K+AFO** or an **HKAFO** into **HK+AFO**. **Compatible with Liber and Indi-Strecker knee joints.**



Weight Classification (body weight)

Measurements

Article No.	AFO/KAFO	Width	Thickness	Length
Separator-S-klein-22	up to 15 kg	16 mm	6,0 + 6,0 mm	80 mm
Separator-S-klein-25-Plus	up to 30 kg	16 mm	6,0 + 6,0 mm	80 mm
Separator-S-mittel-28	up to 40 kg	19 mm	5,5 + 6,0 mm	90 mm
Separator-S-mittel-30	up to 50 kg	19 mm	5,5 + 8,0 mm	90 mm
Separator-S-mittel-30-Plus	up to 70 kg	19 mm	5,5 + 8,0 mm	90 mm
Separator-S-groß-34	up to 110 kg	22 mm	6,0 + 10,0 mm	110 mm



Modular



Detachable




Locked



No Grinding Necessary



Growable



Water Resistant



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com