

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Varo-KA-0-Dynamic-0-N



Features

+ **Force unit** assists knee extension and helps when sitting down and standing up

+ **2-piece adjustable body** to vary flexion stop between 95° and 115°

+ **7 positions** to shift lever arm and increase or decrease torque

Cushioned flexion stop

Wide range of **force options** available

2 different styles of base plate for versatility: retrofit or incorporated




Fig. left


Standard

Fully assembled complete with everything needed for a functional add-on to the Varo-KA-0 knee joint


Available **separately or as combi-product** together with the Varo-KA-0 knee joint

Serial number for traceability

Clearly illustrated instructions for use



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product	Weight Classification (body weight)		Article No.
Varo-KA-0-Dynamic-0-N-28	KAFO	HKAFO/RGO	Dynamic-0-N-28
Varo-KA-0-Dynamic-0-N-30	not yet available		Dynamic-0-N-30
Varo-KA-0-Dynamic-0-N-32	not yet available		Dynamic-0-N-32
Varo-KA-0-Dynamic-0-N-36	up to 50 kg	up to 40 kg	Dynamic-0-N-36
	up to 70 kg	up to 60 kg	

Varo-KA-0 contoured | Varo-KA-0 straight



Locked



Extension




Flexion

Posterior Offset



Modular



Force Options



No Grinding Necessary



Locking Mechanism

Water Resistant



9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com