

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Easy-A and Easy-T Ankle Joint

Features

Waterproof for all situations

- + Active support during **midstance**
- + Active support of the **heel rocker**
- + Resists foot drop during **swing phase**
- + Combine between 2 different **spring forces** or **stop pins**

Noise dampening, quiet operation promotes acceptance

+ Upper section fabricated from **high-strength aluminum or titanium**

Robust, **stainless steel** springs and ball bearings

Bronze movement stops

Range of motion up to **17° dorsiflexion** and **17° plantarflexion** for adjustments to alignment and stops

Exclusive contour for **more stability**

Anchoring points to secure lamination materials

Pre-bent on request

Standard

Fully assembled for clear cost estimates and uncomplicated ordering

Available with **long or short upright**

Upper section **CNC-machined from one piece** of material for increased stability

Machined upright and offset connection ensure **optimal transfer of force**

Precise assembly calibrated to 1/100mm with self-lubricating washers for smooth operation free of play

Maintenance-free stainless steel joint assembly with bronze bushing

Unrivalled **20-month guarantee**

Engraved **serial number** for traceability

Clearly illustrated instructions for use








Fig. left



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)				Joint			
Article No.	Article No.	AFO/KAFO	For More Force		Width	Thickness	Height
Easy-A-19	Easy-T-19	up to 10 kg	up to 10 kg		34 mm	9,0 mm	42 mm
Easy-A-22	Easy-T-22	up to 15 kg	up to 10 kg		39 mm	9,5 mm	46 mm
Easy-A-25	Easy-T-25	up to 25 kg	up to 15 kg		44 mm	11,5 mm	52 mm
Easy-A-28	Easy-T-28	up to 40 kg	up to 30 kg		49 mm	12,0 mm	58 mm

	Modular					Interchangeable Stops			
Adjustable Alignment		Force Options	Range of Motion	Swing Phase	Quick Adjustment		No Grinding Necessary	Noise Dampening	Water Resistant



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com