

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026


Ela-FI-Ex-Post



Features

Changing position of the ball joint
adjusts plantarflexion angle

+ Force unit is used to **influence dorsiflexion**




Standard

Fully assembled complete with everything needed
for a functional add-on to the Ela-FI-Ex ankle joint


Serial number for traceability

Clearly illustrated instructions for use

Fig. left




We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		AFO		Article No.
Ela-FI-Ex-Post-22		up to 15 kg		Post-22
Ela-FI-Ex-Post-25		up to 25 kg		Post-25
Ela-FI-Ex-Post-28		up to 35 kg		Post-28

Adjustable Alignment									Water Resistant
----------------------	---	---	---	--	---	---	---	---	-----------------



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com