

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt


Orthotics 14.0
Valid from July 2026

Shuttle-Flex



Features

- + Tension bands provide a **dampening** function at heel strike and smooth heel rocker. During standing, they provide a **centering**, stabilizing force



Segment attaches to the Shuttle ankle joint stirrup

Standard


Available **separately or as combi-product** together with the Shuttle ankle joint

Precise assembly calibrated to 1/100mm for smooth operation free of play


Serial number for traceability

Clearly illustrated instructions for use

Fig. left









We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		AFO/KAFO	HKAFO/RGO		Article No.
Shuttle-Flex-22		up to 15 kg	up to 15 kg		SH-Flex-22
Shuttle-Flex-25		up to 25 kg	up to 25 kg		SH-Flex-25
Shuttle-Flex-28		up to 35 kg	up to 35 kg		SH-Flex-28
Shuttle-Flex-30		up to 45 kg	up to 45 kg		SH-Flex-30
Shuttle-Flex-32		up to 55 kg	up to 55 kg		SH-Flex-32
Shuttle-Flex-34		up to 70 kg	up to 70 kg		SH-Flex-34





Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com