

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Shuttle-Uno-back



Features

Posterior power unit uses medium to strong forces (plantarflexion resistance) for **reducing** correctable equinus. Weaker forces are used to dampen heel strike for **a smooth heel rocker**

Tension bands provide **a smooth rollover** during stance phase

Segment attaches to the Shuttle ankle joint stirrup




Fig. left


Standard

Available **separately or as combi-product** together with the Shuttle ankle joint

Precise assembly calibrated to 1/100mm for smooth operation free of play

Serial number for traceability


Clearly illustrated instructions for use




We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		AFO/KAFO		Article No.
Shuttle-Uno-back-22		up to 15 kg		SH-Uno-back-22
Shuttle-Uno-back-25		up to 25 kg		SH-Uno-back-25
Shuttle-Uno-back-28		up to 35 kg		SH-Uno-back-28
Shuttle-Uno-back-30		up to 45 kg		SH-Uno-back-30
Shuttle-Uno-back-32		up to 55 kg		SH-Uno-back-32
Shuttle-Uno-back-34		up to 70 kg		SH-Uno-back-34

 Adjustable Alignment	 Modular	 Force Options	 Range of Motion	 Swing Phase	Quick Adjustment	 Interchangeable Stops	 No Grinding Necessary	 Noise Dampening	Water Resistant
---	--	--	--	---	-------------------------	--	--	--	------------------------



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com