

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Shuttle-Uno-front



Features

+ **Anterior power unit** uses interchangeable spring forces (dorsiflexion resistance) **to support** calf musculature

+ **Tension bands** provide **dorsiflexion assist** during swing phase and dampen heel strike for a **smooth heel rocker**

Segment attaches to the Shuttle ankle joint stirrup



Fig. left


Standard

Available **separately or as combi-product** together with the Shuttle ankle joint


Precise assembly calibrated to 1/100mm for smooth operation free of play

Serial number for traceability

Clearly illustrated instructions for use







We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)


Article No. Combi-product		AFO/KAFO		Article No.
Shuttle-Uno-front-22		up to 15 kg		SH-Uno-front-22
Shuttle-Uno-front-25		up to 25 kg		SH-Uno-front-25
Shuttle-Uno-front-28		up to 35 kg		SH-Uno-front-28
Shuttle-Uno-front-30		up to 45 kg		SH-Uno-front-30
Shuttle-Uno-front-32		up to 55 kg		SH-Uno-front-32
Shuttle-Uno-front-34		up to 70 kg		SH-Uno-front-34



ORTHO
SYSTEMS

2000
2026

Years in motion



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com