

**ORTHO**  
SYSTEMS

2000  
2026  
Years in motion



# Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

**Catalogue**  
Excerpt

Orthotics 14.0  
Valid from July 2026

# Indi-Strecker




## Features

Sizes available for the **smallest to large** patients

Changing position of the **ball joint** adjusts flexion angle

+ **Force unit** assists knee extension and helps when sitting down and standing up

Wide range of **force options** available



## Standard


**Fully assembled** complete with everything needed for a functional add-on to the Indi knee joint

Accessory available as **combi-product** with the Indi knee joint


**Serial number** for traceability

**Clearly illustrated** instructions for use


Fig. left





We recommend using our configurator to make your selection.  
The configurator considers, among other factors, weight to height ratios.



#### Weight Classification (body weight)

| Article No. Combi-product |  | KAFO        | HKAFO       |  | Article No.      |
|---------------------------|------------------------------------------------------------------------------------|-------------|-------------|--------------------------------------------------------------------------------------|------------------|
| Indi-Strecker-25-Plus     |                                                                                    | up to 30 kg | up to 25 kg |                                                                                      | Strecker-25-Plus |
| Indi-Strecker-28          |                                                                                    | up to 40 kg | up to 35 kg |                                                                                      | Strecker-28      |
| Indi-Strecker-30          |                                                                                    | up to 50 kg | up to 45 kg |                                                                                      | Strecker-30      |
| Indi-Strecker-30-Plus     |                                                                                    | up to 70 kg | up to 60 kg |                                                                                      | Strecker-30-Plus |
| Indi-Strecker-34          |                                                                                    | up to 85 kg | up to 85 kg |                                                                                      | Strecker-34      |

|                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |  |  |  |  |  |
| Locked                                                                              | Extension 5°                                                                        | Flexion 115°                                                                        | Posterior Offset                                                                    | Modular                                                                               | Force Options                                                                         | No Grinding Necessary                                                                 | Noise Dampening                                                                       | Water Resistant                                                                       |



## Contact

You can reach us from  
by telephone  
per Fax  
per E-Mail  
Address  
Contact persons

9 am to 5 pm (Germany)  
+49 89 904 7575 0  
+49 89 904 7575 5  
[kontakt@ortho-systems-mail.de](mailto:kontakt@ortho-systems-mail.de)  
Merowingerstr. 14, 85551 Kirchheim, Germany  
Jeff Matwichuk and Thomas Böckh

[www.ortho-systems-inmotion.com](http://www.ortho-systems-inmotion.com)