

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Klassik-Konso Knee Joint



Features

Sizes available for the **smallest to large** patients

Joint flexes to **130°**

Fabricated from **high-strength aluminum**

+ Adjustable flexion angle from **0° to 40°** in 5° intervals

Available with a contoured or straight lower bar

Standard

Fully assembled for clear cost estimates and uncomplicated ordering

Integrated upright bars ensure **optimal transfer of force**

Upper section **CNC-machined from one piece** of material for increased stability

Precise assembly calibrated to 1/100mm with self-lubricating washers for smooth operation free of play

Maintenance-free stainless steel joint assembly with bronze bushing

Unrivalled **20-month guarantee**

Engraved **serial number** for traceability

Clearly illustrated instructions for use





Fig. left




Contact


You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons


9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)			Joint			
Article No.	KAFO	HKAFO/RGO		Width	Thickness	Thickness at Drop Lock
Klassik-Konso-22	up to 10 kg	up to 10 kg		22 mm	10,5 mm	15,0 mm
Klassik-Konso-25	up to 20 kg	up to 5 kg		25 mm	13,0 mm	16,5 mm
Klassik-Konso-28	up to 30 kg	up to 20 kg		28 mm	15,0 mm	19,5 mm
Klassik-Konso-30	up to 35 kg	up to 30 kg		30 mm	16,0 mm	20,5 mm
Klassik-Konso-32	up to 50 kg	up to 40 kg		32 mm	17,0 mm	21,0 mm
Klassik-Konso-36	up to 70 kg	up to 55 kg		34 mm	17,0 mm	21,0 mm
Klassik-Konso-40	up to 85 kg	up to 65 kg		40 mm	21,0 mm	26,0 mm
Klassik-Konso-40-Plus	up to 85 kg	up to 85 kg	40 mm	22,0 mm	28,0 mm	

Klassik-Konso contoured | Klassik-Konso straight

								
Locked	Extension	Flexion	Posterior Offset	Modular	Force Options	No Grinding Necessary	Locking Mechanism	Water Resistant