

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Klassik-Konso Slide



Features

+ Wear-free, **magnetic** safety release

+ Ergonomic **plastic grip**

Pre-select system opens only when knee is in full extension and drop lock is relieved

Cable with or without housing




Fig. left


Standard

Fully assembled for clear cost estimates and uncomplicated ordering


Accessory available as **combi-product** with the Klassik-Konso knee joint

Serial number for traceability


Clearly illustrated instructions for use



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		KAFO		Article No.
Klassik-Konso-Slide-30		up to 35 kg		Slide-30
Klassik-Konso-Slide-32		up to 50 kg		Slide-32
Klassik-Konso-Slide-36		up to 70 kg		Slide-36
Klassik-Konso-Slide-40		up to 85 kg		Slide-40
Klassik-Konso-Slide-40-Plus		up to 85 kg		Slide-40-Plus
				Only available as replacement item



Locked



Extension



Flexion

Posterior Offset

Modular

Force Options



No Grinding Necessary



Locking Mechanism

Water Resistant

www.ortho-systems-inmotion.com