

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Konso-0-Move-0-N



Features

- Integrated ball plunger** enables temporary unlocking of the joint
- + Force unit** assists knee extension and helps when sitting down and standing up
- + 2-piece adjustable body** to vary flexion stop between **95°** and **115°**
- + 7 positions** to shift lever arm and increase or decrease torque
- Cushioned** flexion stop
- Wide range of **force options** available
- 2 different styles** of base plate for versatility: retrofit or incorporated




Fig. left


Standard

Fully assembled complete with everything needed for a functional add-on to the Konso knee joint

Available **separately** or as **combi-product** together with a Konso knee joint

Serial number for traceability


Clearly illustrated instructions for use




We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		KAFO	HKAFO/RGO		Article No.
Konso-0-Move-0-N-25		up to 20 kg	up to 15 kg		Move-0-N-25
Konso-0-Move-0-N-28		up to 30 kg	up to 20 kg		Move-0-N-28
Konso-0-Move-0-N-30		up to 35 kg	up to 30 kg		Move-0-N-30
Konso-0-Move-0-N-32		up to 50 kg	up to 40 kg		Move-0-N-32
Konso-0-Move-0-N-36		up to 70 kg	up to 60 kg		Move-0-N-36
Konso-0-Move-0-N-40		up to 85 kg	up to 85 kg		Move-0-N-40

Konso-0 contoured | Konso-0 straight

 Locked	 Extension	 Flexion	 Posterior Offset	 Modular	 Force Options	 No Grinding Necessary	 Locking Mechanism	 Water Resistant
---	--	--	---	--	--	--	--	--



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com