

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Konso-5-Move-5-S



Features




Fig. left


Standard

Fully assembled complete with everything needed for a functional add-on to the Konso knee joint


Available **separately** or as **combi-product** together with a Konso knee joint

Serial number for traceability



Clearly illustrated instructions for use




We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.




Weight Classification (body weight)

Article No. Combi-product			KAFO	HKAFO / RGO		
Konso-5-Move-5-S-25			up to 20 kg	up to 15 kg		Move-5-S-25
Konso-5-Move-5-S-28			up to 30 kg	up to 20 kg		Move-5-S-28
Konso-5-Move-5-S-30			up to 35 kg	up to 30 kg		Move-5-S-30
Konso-5-Move-5-S-32			up to 50 kg	up to 40 kg		Move-5-S-32
Konso-5-Move-5-S-36			up to 70 kg	up to 60 kg		Move-5-S-36
Konso-5-Move-5-S-40			up to 85 kg	up to 85 kg		Move-5-S-40

Konso-5 contoured | Konso-5 straight

 Locked	 Extension	 Flexion	Posterior Offset	 Modular	 Force Options	 No Grinding Necessary	 Locking Mechanism	Water Resistant
---	--	--	------------------	--	--	--	--	-----------------



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com