

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics


Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026


Assembly Tools

Precision tools for secure, accurate, and efficient assembly of our orthotic joints




Montage-A-EH-GV

Assembly aid for all Blocker, Konso, Varo, Klassik-Konso, Liber, Liber-IN, Indi, Shuttle, Ela, Naveta and Jumper joints





Montage-K-GV-kurz

Short pin spanner for all Fokus and Modos hip joints, and size -22 Jumper, Naveta, Ela-FI-Ex and Ela-Dor ankle joints




Montage-EH-GV-(Size)

Assembly aid for Fokus and Modos 3D and 2D hip joints




Montage-G-GV-lang

Long pin spanner for all Fokus and Modos hip joints, and size-22 Jumper, Naveta, Ela-Fl-Ex and Ela-Dor ankle joints




Montage-H-Rezi-WM

Pin spanner for RGO axle and rocker bar connections



Montage-B-Rezi-GV-WM

Hex nut driver for inner rocker bar connection





Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com