

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics

Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Donning Aid

The Donning Aid Rudi can be used with all ORTHO-SYSTEMS extension-assist cylinders. It clips over the cylinder to temporarily fix the position between 90° and 150° (30°). Keeping the knee in a set position can make it easier to put on the orthosis. Additionally, if the cylinder has ball joints on both ends, it can be completely removed together with the donning aid.



Donning Aid Rudi

Liber-Extensor-S fixed at 95° (85°)


Compatible with:

Extensor-S, Extensor-N
Extensor-IN-S, Extensor-IN-N
Strecker
Dynamic-0-S, Dynamic-0-N
Move-0/5, Move-0/5-N

Indi-Strecker with cylinder removed

Compatible with:

Extensor-S
Extensor-IN-S
Strecker
Dynamic-0-S
Move-0/5-S



for Extension Assist

Article No.	for Cylinder Diameter	Cylinder Travel
Anziehhilfe-Rudi-12-50-klein	12 mm	40 and 50 mm
Anziehhilfe-Rudi-15-60-mittel	15 mm	50 and 60 mm
Anziehhilfe-Rudi-18-60-groß	18 mm	60 mm

► See pages 76 - 79 for applications



Donning Aid Rudi

Extensor-	Dynamic-	Move-
22, 25-Plus	-	25
28, 30, 30-Plus	28, 30, 32, 36	28, 30, 32, 36
34	-	40



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com